

## Biopsy Forceps $\varnothing$ 2.3 mm for working channels $\geq \varnothing$ 2.8 mm

stainless steel spiral, smoothed at distal end:



length	oval fenestrated cups	
	without spike	with spike
160 cm	UK 02 22V-B	UK 02 22VD-B
180 cm	UK 02 22V-C	UK 02 22VD-C
230 cm	UK 02 22V-G	UK 02 22VD-G



**These forceps are construction-conditioned non-repairable.**

## HF connection fitting to active cord:



### Hot Biopsy Forceps $\varnothing$ 2.3 mm, for working channels $\geq \varnothing$ 2.8 mm

Teflon® coated stainless steel spiral:



length	oval fenestrated cups
	without spike
180 cm	OK 61 22-C
230 cm	OK 61 22-G



### Hot Biopsy Forceps $\varnothing$ 2.5 mm, for working channels $\geq \varnothing$ 2.8 mm

Teflon® coated stainless steel spiral:



length	oval fenestrated cups
	without spike
180 cm	OK 61 25-C
230 cm	OK 61 25-G



## recommended equipment for Hot Biopsy

active cord:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument possesses a flushing port.



This instrument is designed to be used with HF current.

ENDO-FLEX GmbH  
Alte Hünxer Str. 113 - 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: export@endo-flex.de  
www.endo-flex.de

# BIOPSY FORCEPS

re-usable



 **ENDO-FLEX**®

## Biopsy Forceps $\varnothing$ 2.3 mm for working channels $\geq \varnothing$ 2.8 mm

stainless steel spiral, smoothed at distal end:




length	round fenestrated cups		oval fenestrated cups	
	without spike	with spike	without spike	with spike
120 cm	K 01 22V-A	K 01 22VD-A	K 02 22V-A	K 02 22VD-A
160 cm	K 01 22V-B	K 01 22VD-B	K 02 22V-B	K 02 22VD-B
180 cm	K 01 22V-C	K 01 22VD-C	K 02 22V-C	K 02 22VD-C
230 cm	K 01 22V-G	K 01 22VD-G	K 02 22V-G	K 02 22VD-G

stainless steel spiral, smoothed at distal end:



length	oval toothed cups		oval cups rat tooth
	without spike	with spike	without spike
120 cm	K 04 22V-A	K 04 22VD-A	K 07 22V-A
160 cm	K 04 22V-B	K 04 22VD-B	K 07 22V-B
180 cm	K 04 22V-C	K 04 22VD-C	K 07 22V-C
230 cm	K 04 22V-G	K 04 22VD-G	K 07 22V-G

overall stainless smoothed steel spiral:



length	oval fenestrated cups	
	without spike	with spike
120 cm	KF 02 22-A	KF 02 22D-A
160 cm	KF 02 22-B	KF 02 22D-B
180 cm	KF 02 22-C	KF 02 22D-C
230 cm	KF 02 22-G	KF 02 22D-G

stainless steel spiral, smoothed at distal end, cups tiltable in both directions:



length	oval fenestrated cups
	without spike
120 cm	K 02 22V-AK
160 cm	K 02 22V-BK
180 cm	K 02 22V-CK
230 cm	K 02 22V-GK



## Biopsy Forceps $\varnothing$ 3.4 mm for working channels $\geq \varnothing$ 3.7 mm

stainless steel spiral, smoothed at distal end:



length	round fenestrated cups	
	without spike	with spike
120 cm	K 01 34V-A	K 01 34VD-A
160 cm	K 01 34V-B	K 01 34VD-B
180 cm	K 01 34V-C	K 01 34VD-C
230 cm	K 01 34V-G	K 01 34VD-G



## Biopsy Forceps $\varnothing$ 2.5 mm for working channels $\geq \varnothing$ 2.8 mm

stainless steel spiral, smoothed at distal end:



length	oval fenestrated cups	
	without spike	with spike
120 cm	K 02 25V-A	K 02 25VD-A
160 cm	K 02 25V-B	K 02 25VD-B
180 cm	K 02 25V-C	K 02 25VD-C
230 cm	K 02 25V-G	K 02 25VD-G



## Biopsy Forceps $\varnothing$ 2.1 mm for working channels $\geq \varnothing$ 2.4 mm

stainless steel spiral, smoothed at distal end:



length	short oval fenestrated cups	
	without spike	with spike
120 cm	K 02 21V-A	K 02 21VD-A
160 cm	K 02 21V-B	K 02 21VD-B
180 cm	K 02 21V-C	K 02 21VD-C
230 cm	K 02 21V-G	K 02 21VD-G



## Biopsy Forceps $\varnothing$ 1.8 mm for working channels $\geq \varnothing$ 2.0 mm

stainless steel spiral, smoothed at distal end:



length	oval fenestrated cups		short oval fenestrated cups	
	without spike	with spike	without spike	with spike
120 cm	K 16 18V-A	K 16 18VD-A	K 17 18V-A	K 17 18VD-A
160 cm	K 16 18V-B	K 16 18VD-B	K 17 18V-B	K 17 18VD-B
180 cm	K 16 18V-C	K 16 18VD-C	K 17 18V-C	K 17 18VD-C
230 cm	K 16 18V-G	K 16 18VD-G	K 17 18V-G	K 17 18VD-G

stainless steel spiral, smoothed at distal end:



length	oval toothed cups
	without spike
120 cm	K 04 18V-A
160 cm	K 04 18V-B
180 cm	K 04 18V-C
230 cm	K 04 18V-G



## Mini Biopsy Forceps for working channels $\geq \varnothing$ 1.2 mm

overall smoothed stainless steel spiral:



length	oval cups	
	$\varnothing$ 0.8 mm	$\varnothing$ 1.0 mm
120 cm		KF 02 10-A
160 cm		KF 02 10-B
180 cm		KF 02 10-C
230 cm	KF 02 08-G	KF 02 10-G

