

MUCOSECTOMY SNARES

disposable



 **ENDO-FLEX[®]**


HF connection fitting to active cord:



disposable Mucosectomy Snare

Ø 2.3 mm, plastic handle, for working channels \geq Ø 2.8 mm


oval loop



length	package	25 mm
230 cm	5 pieces	OE 34 2246-G

recommended equipment

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

Legend:



This instrument is designed and warranted for single use only.



This instrument is designed to be used with HF current.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
 Alte Hünxer Str. 113 - 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de