

# POLYPECTOMY SNARES

disposable



 **ENDO-FLEX**<sup>®</sup>


# HF connection fitting to active cord:

rotatable disposable Polypectomy Snare

Ø 2.3 mm, plastic handle, for working channels ≥ Ø 2.8 mm



oval loop braided wire



length	package	15 mm	25 mm	35 mm
180 cm	5 pieces	NOE 34 2214-C	NOE 34 2216-C	NOE 34 2217-C
230 cm	5 pieces	NOE 34 2214-G	NOE 34 2216-G	NOE 34 2217-G


disposable Polypectomy Snare Ø 2.3 mm, plastic handle, for working channels ≥ Ø 2.8 mm

oval loop braided wire



length	package	20 mm	25 mm	35 mm
180 cm	5 pieces	OE 34 2215-C	OE 34 2216-C	OE 34 2217-C
230 cm	5 pieces	OE 34 2215-G	OE 34 2216-G	OE 34 2217-G

oval loop monofilament wire



length	package	25 mm
180 cm	5 pieces	OE 34 2216M-C
230 cm	5 pieces	OE 34 2216M-G





## recommended equipment

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

### Legend:

-  This instrument is designed and warranted for single use only.
-  This instrument is designed to be used with HF current.
-  This instrument is rotatable.
-  This instrument is sterile packed separately.

ENDO-FLEX GmbH  
 Alte Hünxer Str. 113 - 115  
 46562 Voerde / Germany  
 Phone: +49 281 - 9 44 00 - 0  
 Fax: +49 281 - 9 44 00 - 11  
 e-mail: export@endo-flex.de  
 www.endo-flex.de