

MUCOSECTOMY SNARES

re-usable


AUTOCLEAVABLE
134°C
273°F



 **ENDO-FLEX**[®]

rotatable Mucosectomy Snare \varnothing 2.3 mm,
without handle, for working channels $\geq \varnothing$ 2,8 mm

oval loop



monofilament wire

length	20 mm	25 mm	35 mm
180 cm	D 34 22151M-C	D 34 22161M-C	D 34 22171M-C
230 cm	D 34 22151M-G	D 34 22161M-G	D 34 22171M-G

recommended equipment

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

universal handle:

Order No.:	62 0300
------------	---------



Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.



This instrument is designed to be used with HF current.



This instrument is rotatable.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de