


Polypectomy Snare \varnothing 2.3 mm, incl. handle,
for working channels $\geq \varnothing$ 2.8 mm

asymmetrical loop



length	40 mm		60 mm	
	standard tip	insulated tip	standard tip	insulated tip
180 cm	OK 65 2240-C	OK 65 22401-C	OK 65 2260-C	OK 65 22601-C
230 cm	OK 65 2240-G	OK 65 22401-G	OK 65 2260-G	OK 65 22601-G

Polypectomy Snare \varnothing 2.3 mm, fixed handle,
for working channels $\geq \varnothing$ 2.8 mm \varnothing



length	35 mm	20 - 35 mm
	standard tip	standard tip
230 cm	OK 34 2237-G	OK 34 2267-G

recommended equipment

active cords:




HF generators	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500





universal handle:

(for snares, which are supplied without handle)



Order No.:	62 0300
------------	---------



Legend:

-  This instrument is autoclavable by 134°C / 273°F.
-  This instrument is detachable for cleaning.
-  This instrument is designed to be used with HF current.
-  This Instrument possesses a flushing port.

POLYPECTOMY SNARES
re-usable



© ENDO-FLEX - Technical specifications are subject to change

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de




HF connection fitting to active cord:

Polypectomy Snare \varnothing 2.3 mm, fixed handle, for working channels $\geq \varnothing$ 2.8 mm




asymmetrical loop




length	15 mm		25 mm	
	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 30 22-C	OK 32 22-C	OK 31 22-C	OK 33 22-C
230 cm	OK 30 22-G	OK 32 22-G	OK 31 22-G	OK 33 22-G

oval loop



length	20 mm		25 mm		35 mm	
	insulated tip	standard tip	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 34 2205-C	OK 34 2215-C	OK 34 2206-C	OK 34 2216-C	OK 34 2207-C	OK 34 2217-C
230 cm	OK 34 2205-G	OK 34 2215-G	OK 34 2206-G	OK 34 2216-G	OK 34 2207-G	OK 34 2217-G

hexagonal loop



length	20 mm		25 mm		35 mm	
	insulated tip	standard tip	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 38 2205-C	OK 37 2205-C	OK 38 2206-C	OK 37 2206-C	OK 38 2207-C	OK 37 2207-C
230 cm	OK 38 2205-G	OK 37 2205-G	OK 38 2206-G	OK 37 2206-G	OK 38 2207-G	OK 37 2207-G

Polypectomy Snare \varnothing 2.3 mm, fixed handle, for working channels $\geq \varnothing$ 2.8 mm


oval loop monofilament wire



length	20 mm	25 mm	35 mm
	standard tip	standard tip	standard tip
180 cm	OK 34 2215M-C	OK 34 2216M-C	OK 34 2217M-C
230 cm	OK 34 2215M-G	OK 34 2216M-G	OK 34 2217M-G

Polypectomy Snare \varnothing 2.3 mm, without handle, for working channels $\geq \varnothing$ 2.8 mm

asymmetrical loop



length	15 mm		25 mm		15 mm		25 mm	
	insulated tip		insulated tip		standard tip		standard tip	
	standard	rotatable	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	30 22-C	D 30 22-C	31 22-C	D 31 22-C	32 22-C	D 32 22-C	33 22-C	D 33 22-C
230 cm	30 22-G	D 30 22-G	31 22-G	D 31 22-G	32 22-G	D 32 22-G	33 22-G	D 33 22-G

insulated tip oval loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2205-C	D 34 2205-C	34 2206-C	D 34 2206-C	34 2207-C	D 34 2207-C
230 cm	34 2205-G	D 34 2205-G	34 2206-G	D 34 2206-G	34 2207-G	D 34 2207-G

standard tip oval loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2215-C	D 34 2215-C	34 2216-C	D 34 2216-C	34 2217-C	D 34 2217-C
230 cm	34 2215-G	D 34 2215-G	34 2216-G	D 34 2216-G	34 2217-G	D 34 2217-G

insulated tip hexagonal loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	38 2205-C	D 38 2205-C	38 2206-C	D 38 2206-C	38 2207-C	D 38 2207-C
230 cm	38 2205-G	D 38 2205-G	38 2206-G	D 38 2206-G	38 2207-G	D 38 2207-G

necessary devices: universal handle 62 0300



Polypectomy Snare \varnothing 2.3 mm, without handle, for working channels $\geq \varnothing$ 2.8 mm

standard tip hexagonal loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	37 2205-C	D 37 2205-C	37 2206-C	D 37 2206-C	37 2207-C	D 37 2207-C
230 cm	37 2205-G	D 37 2205-G	37 2206-G	D 37 2206-G	37 2207-G	D 37 2207-G

standard tip oval loop



monofilament wire




length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2215M-C	D 34 2215M-C	34 2216M-C	D 34 2216M-C	34 2217M-C	D 34 2217M-C
230 cm	34 2215M-G	D 34 2215M-G	34 2216M-G	D 34 2216M-G	34 2217M-G	D 34 2217M-G

Polypectomy Snare \varnothing 1.8 mm, without handle, for working channels $\geq \varnothing$ 2.0 mm

oval loop



length	20 mm	25 mm	20 mm	25 mm
	insulated tip		standard tip	
180 cm	34 1805-C	34 1806-C	34 1815-C	34 1816-C
230 cm	34 1805-G	34 1806-G	34 1815-G	34 1816-G



length	15 mm	25 mm	20 mm	25 mm	20 mm	25 mm
	asymmetrical loop, standard tip		hexagonal loop, insulated tip		hexagonal loop, standard tip	
180 cm	32 18-C	33 18-C	38 1805-C	38 1806-C	37 1805-C	37 1806-C
230 cm	32 18-G	33 18-G	38 1805-G	38 1806-G	37 1805-G	37 1806-G

necessary devices: universal handle 62 0300

