


FIBRIN APPLICATION NEEDLES

disposable



 **ENDO-FLEX**[®]

disposable double lumen probes **Ø 2.6 mm** for Fibrin application,
for working channels $\geq \text{Ø } 2.8 \text{ mm}$

needle Ø 0.9 mm,
needle length: 6 mm



Particularly designed to be used along with
Beriplast® HS Combi-Set.

- binding of both components in tissue
- simultaneously injection of both components



incl. needle safety device

- colour marked LL connectors:

RED = Trombin

BLUE = Fibrinogen

length	Order No.
160 cm	ET 27 2620-B
180 cm	ET 27 2620-C
230 cm	ET 27 2620-G

coaxial needle Ø 1.2 mm,
needle length: 6 mm



- binding of both components inside of tissue
- simultaneously injection of both components



incl. needle safety device

- colour marked LL connectors:

RED = Trombin

BLUE = Fibrinogen

length	Order No.
160 cm	ET 27 2630-B
180 cm	ET 27 2630-C
230 cm	ET 27 2630-G

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de