

SCLEROTHERAPY NEEDLES

disposable

A blurred background image showing a person in a light blue medical uniform standing next to a medical cart or equipment in a clinical setting.


 **ENDO-FLEX**[®]




disposable Sclerotherapy Needles \varnothing 2.3 mm, for working channels $\geq \varnothing$ 2.8 mm

needle \varnothing 0.7 mm

needle \varnothing 0.9 mm





incl. needle safety device 




length	package	needle \varnothing 0.7 mm		needle \varnothing 0.9 mm	
		needle length		needle length	
		4 mm	6 mm	4 mm	6 mm
160 cm	5 pieces		ET 25 22-B6		ET 26 22-B6
180 cm	5 pieces	ET 25 22-C4	ET 25 22-C6		ET 26 22-C6
230 cm	5 pieces	ET 25 22-G4	ET 25 22-G6		ET 26 22-G6

disposable Sclerotherapy Needles \varnothing 1.8 mm, for working channels $\geq \varnothing$ 2.0 mm

needle 0.6 mm \varnothing



incl. needle safety device 

length	package	needle \varnothing 0,6 mm	
		needle length	
		4 mm	6 mm
180 cm	5 pieces	ET 25 18-C4	ET 25 18-C6

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de