

SCLEROTHERAPY NEEDLES

re-usable

AUTOCLEAVABLE
134°C
273°F



 **ENDO-FLEX**[®]

Sclerotherapy Needles \varnothing 2.3 mm, TEFLON[®]-Tube, for working channels $\geq \varnothing$ 2.8 mm


needle \varnothing 0.7 mm



length	standard needle length		flushable needle length	
	4 mm	6 mm	4 mm	6 mm
160 cm	T 25 22-B4	T 25 22-B6	T 25 22S-B4	T 25 22S-B6
180 cm	T 25 22-C4	T 25 22-C6	T 25 22S-C4	T 25 22S-C6
230 cm	T 25 22-G4	T 25 22-G6	T 25 22S-G4	T 25 22S-G6

Sclerotherapy Needles \varnothing 1.8 mm, TEFLON[®]-Tube, for working channels $\geq \varnothing$ 2.0 mm

needle \varnothing 0.6 mm



length	standard needle length	
	4 mm	6 mm
160 cm	T 25 18-B4	T 25 18-B6
180 cm	T 25 18-C4	T 25 18-C6
230 cm	T 25 18-G4	

Legend:



This instrument is autoclavable by 134°C / 273°F.

This instrument is detachable for cleaning.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de