

CLEANING BRUSHES

disposable



 **ENDO-FLEX[®]**

disposable Cleaning Brushes \varnothing 1.8 mm for cleaning of endoscope working channels

double sided,
length of brushes 20 mm



length	package	\varnothing 1.8 mm brush \varnothing 5.0 mm	\varnothing 1.8 mm brush \varnothing 5.0 mm
		for working channels $\geq \varnothing$ 2.0 mm	for working channels $\geq \varnothing$ 2.8 mm
120 cm	10 pieces	EBS 41 181-A	
180 cm	10 pieces		EBS 41 221-C
230 cm	10 pieces		EBS 41 221-G

disposable Cleaning Brushes for cleaning of endoscope valves

double sided



small brush \varnothing	large brush \varnothing	Order No.
4.0 mm	12.0 mm	ZB0040

Legend:



This instrument is designed and warranted for single use only.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de