

# CLEANING BRUSHES

re-usable

AUTOCLEANABLE  
134°C  
273°F



 **ENDO-FLEX<sup>®</sup>**

# Brushes $\varnothing$ 1.8 and 2.2 mm for cleaning of endoscope working channels

single sided, incl. handle



length	spiral $\varnothing$ 1.8 mm brush $\varnothing$ 3.0 mm	spiral $\varnothing$ 2.2 mm brush $\varnothing$ 5.0 mm	spiral $\varnothing$ 2.2 mm brush $\varnothing$ 7.0 mm
	for working channels $\geq \varnothing$ 2.0 mm	for working channels $\geq \varnothing$ 2.8 mm	for working channels $\geq \varnothing$ 3.2 mm
120 cm	41 181-A	41 221-A	41 401-A
160 cm	41 181-B	41 221-B	41 401-B
180 cm	41 181-C	41 221-C	41 401-C
230 cm		41 221-G	41 401-G

double sided



length	spiral $\varnothing$ 1.8 mm brush $\varnothing$ 3.0 mm	spiral $\varnothing$ 2.2 mm brush $\varnothing$ 5.0 mm	spiral $\varnothing$ 2.2 mm brush $\varnothing$ 7.0 mm
	for working channels $\geq \varnothing$ 2.0 mm	for working channels $\geq \varnothing$ 2.8 mm	for working channels $\geq \varnothing$ 3.2 mm
120 cm	BS 41 181-A	BS 41 221-A	BS 41 401-A
160 cm	BS 41 181-B	BS 41 221-B	BS 41 401-B
180 cm	BS 41 181-C	BS 41 221-C	BS 41 401-C
230 cm		BS 41 221-G	BS 41 401-G

Legend:



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH  
Alte Hünxer Str. 113 - 115  
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)