



E.R.C.P. Catheters  $\varnothing$  1.8 mm for working channels  $\geq \varnothing$  2.0 mm,  
not to be used with guide wire


needle tip 2 mm




type	3 black markers	coloured marker	connector
with mandrin	49 18	49 181	

E.R.C.P. Catheters  $\varnothing$  1.6 mm for working channels  $\geq \varnothing$  2.0 mm,  
accepts guide wire .025 inch

tapered tip



execution	3 black markers	coloured markers	connector
with mandrin	44 16	44 161	

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400



© ENDO-FLEX - Technical specifications are subject to change

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is sterile packed separately.

ENDO-FLEX GmbH  
Alte Hünxer Str. 113 - 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: export@endo-flex.de  
www.endo-flex.de


**E.R.C.P. CATHETERS**  
re-usable



 **ENDO-FLEX**<sup>®</sup>


E.R.C.P. Catheters  $\varnothing$  1.8 mm,  
accepts guide wire .035 inch

standard tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	43 18	43 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	43 22-18FD	43 221-18FD	


tapered tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	44 18	44 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	44 22-18FD	44 221-18FD	


E.R.C.P. Catheters  $\varnothing$  1.8 mm,  
accepts guide wire .035 inch

conical  
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	45 18	45 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	45 22-18FD	45 221-18FD	

tapered  
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	46 18	46 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	46 22-18FD	46 221-18FD	


E.R.C.P. Catheters  $\varnothing$  1.8 mm,  
accepts guide wire .035 inch

round  
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	47 18	47 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	47 22-18FD	47 221-18FD	

hour glass  
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels $\geq \varnothing$ 2.0 mm	48 18	48 181	
with mandrin / proximal $\varnothing$ 2.45 mm for working channels $\geq \varnothing$ 2.8 mm	48 22-18FD	48 221-18FD	