

# COAGULATION PROBES

re-usable



**ENDO-FLEX**<sup>®</sup>

## HF connector fitting to active cord:

### Coagulation Probe $\varnothing$ 1.8 and 2.3 mm



Tip with 1 opening,  
handle with flushing port



length	$\varnothing$ 1.8 mm for working channels $\geq \varnothing$ 2.0 mm	$\varnothing$ 2.3 mm for working channels $\geq \varnothing$ 2.8 mm
120 cm	901 80 18-A	901 80 22-A
160 cm	901 80 18-B	901 80 22-B
180 cm	901 80 18-C	901 80 22-C
230 cm	901 80 18-G	901 80 22-G

## recommended equipment

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

### Legend:



This instrument is autoclavable by 134°C / 273°F.



This Instrument possesses a flushing port.



This instrument is designed to be used with HF current.

ENDO-FLEX GmbH  
 Alte Hünxer Str. 113 - 115  
 46562 Voerde / Germany  
 Phone: +49 281 - 9 44 00 - 0  
 Fax: +49 281 - 9 44 00 - 11  
 e-mail: export@endo-flex.de  
 www.endo-flex.de