

CYSTO GASTRO SET

re-usable


AUTOCCLAVABLE
134°C
273°F



 **ENDO-FLEX**®

Cysto Gastro Set consisting of needle and catheter

HF-connection fitting
to active cord:



length	Fr. 8.5	Fr. 10.0
210 cm	XS99900901341	XS99900901340

equipment

active cord:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is designed to be used with HF current.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de