



disposable Needle-Knife Sphincterotomes \varnothing 1.8 mm, fixed handle,
for working channels $\geq \varnothing$ 2.0 mm


HF connection fitting to active cord:



pointed type,
stop at 3 mm

 single lumen

OE 110 18N3




blunt type,
stop at 3 mm

 single lumen

OE 110 18SP3

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

Legend:



This instrument is designed and warranted for single use only.



This instrument is designed to be used with HF current.



This instrument is sterile packed separately.




This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de

SPHINCTEROTOMES


disposable



 ENDO-FLEX®


disposable **non-guide-wire** Sphincterotomes **Ø 1.8 mm**, fixed handle, for working channels $\geq \text{Ø } 2.0 \text{ mm}$

HF connection fitting to active cord:




tapered tip

	 single lumen	
cutting length	20 mm	30 mm
	OE 105 1821	OE 105 1831



precut type

	 single lumen		 single lumen	
cutting length	20 mm	30 mm	20 mm	30 mm
	OE 103 1821	OE 103 1831	OE 105 1821PC	OE 105 1831PC



precut type with tapered tip

disposable **guide wire** Sphincterotomes, **Ø 2.2 mm** distal **Ø 1.8 mm**, fixed handle, for working channels $\geq 2.8 \text{ mm } \text{Ø}$

HF connection fitting to active cord:





standard tip

	 double lumen		
cutting length	20 mm	25 mm	30 mm
	OE 100 2220DL	OE 100 2225DL	OE 100 2230DL



precut type

	 double lumen			 double lumen		
cutting length	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
	OE 102 2220DL	OE 102 2225DL	OE 102 2230DL	OE 104 2220DL	OE 104 2225DL	OE 104 2230DL



tapered tip

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400



disposable **guide wire** Sphincterotomes, **Ø 2.2 mm** distal **Ø 1.8 mm**, fixed handle, for working channels $\geq 2.8 \text{ mm } \text{Ø}$

HF connection fitting to active cord:



precut type

	 triple lumen		
cutting length	20 mm	25 mm	30 mm
	OE 102 2220TRL	OE 102 2225TRL	OE 102 2230TRL



tapered tip

	 triple lumen		
cutting length	20 mm	25 mm	30 mm
	OE 104 2220TRL	OE 104 2225TRL	OE 104 2230TRL

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400

