



Sphincterotomes  $\varnothing$  1.8 mm, fixed handle,  
for working channels  $\geq \varnothing$  2.0 mm

HF connection fitting to active cord:





	<b>single lumen</b>	<b>single lumen</b>
standard	OK 110 18SP	OK 110 18N



	<b>single lumen</b>
standard	OK 112 18

Sphincterotomes  $\varnothing$  2.2 mm, fixed handle,  
for working channels  $\geq \varnothing$  2.8 mm

HF connection fitting to active cord:



	<b>single lumen</b>
standard	OK 112 22
accepts guide wire .035 inch	OK 112 22FD

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400

active cords:

HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500


Legend:

- This instrument is autoclavable by 134°C / 273°F.
- This instrument is designed to be used with HF current.
- This Instrument possesses a flushing port.
- This instrument is sterile packed separately.

ENDO-FLEX GmbH  
Alte Hünxer Str. 113 - 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: export@endo-flex.de  
www.endo-flex.de

# SPHINCTEROTOMES

re-usable



© ENDO-FLEX - technical specifications are subject to change

non-guide-wire Sphincterotomes  $\varnothing$  1.8 mm, fixed handle, for working channels  $\geq$  2.0 mm  $\varnothing$

HF connection fitting to active cord:





	single lumen		single lumen	
cutting length	20 mm	30 mm	20 mm	30 mm
standard	OK 101 1820	OK 101 1830	OK 103 1820	OK 103 1830
cutting direction 11 o'clock	OK 101 1822	OK 101 1832	OK 103 1822	OK 103 1832



	single lumen		single lumen	
cutting length	20 mm	30 mm	20 mm	30 mm
standard	OK 105 1820	OK 105 1830	OK 105 1820PC	OK 105 1830PC
cutting direction 11 o'clock	OK 105 1822	OK 105 1832	OK 105 1822PC	OK 105 1832PC

single lumen guide wire Sphincterotomes,  $\varnothing$  2.2 mm distal  $\varnothing$  1.8 mm, fixed handle, for working channels  $\geq$   $\varnothing$  2.8 mm

HF connection fitting to active cord:



	single lumen		single lumen	
cutting length	20 mm	30 mm	20 mm	30 mm
standard	OK 103 1820GW	OK 103 1830GW	OK 105 1820GW	OK 105 1830GW
cutting direction 11 o'clock	OK 103 1820GW11	OK 103 1830GW11	OK 105 1820GW11	OK 105 1830GW11



	double lumen			double lumen		
cutting length	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
standard	OK 100 2220DL	OK 100 2225DL	OK 100 2230DL	OK 103 2220DL	OK 103 2225DL	OK 103 2230DL

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400



double lumen guide wire Sphincterotomes,  $\varnothing$  2.2 mm distal  $\varnothing$  1.8 mm, fixed handle, for working channels  $\geq$   $\varnothing$  2.8 mm

HF connection fitting to active cord:



	double lumen			double lumen
cutting length	20 mm	25 mm	30 mm	—
standard	OK 105 2220DL	OK 105 2225DL	OK 105 2230DL	OK 110 18N-DL

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400

