


BILIARY DILATION CATHETERS

re-usable



Biliary Dilation Catheters

- catheter length 180 cm
- radiopaque teflon tube
- distal radiopaque markers



catheter tip	Fr. 3.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 7.0	Fr. 7.0
catheter Ø	Fr. 5.0	Fr. 7.0	Fr. 8.5	Fr. 9.0	Fr. 10.0	Fr. 10.0	Fr. 12.0
accepts guide wire	.025 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch
Order No.	200 00053	200 00075	200 00085	200 00095	200 00105	200 00107	200 00127

recommended equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400



Legend:



This Instrument possesses a flushing port.



This instrument is sterile packed separately.



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de