


P T C D
disposable



 **ENDO-FLEX[®]**

disposable PTC-Set Fr. 7.0 / consisting of 4 parts

PTC needle



ring-drainage-catheter



guide wire



PTC exchange guide wire



Order No.	name
280 170	PTCD-Set Fr. 7.0 / consisting of 4 parts

single components:

280 171	PTC needle 19 G incl. catheter Fr. 5.0 / length 30 cm
280 172	teflon coated stainless steel guide wire .035 inch / length 100 cm
280 173	PTC exchanging guide wire .035 inch / 7.5 cm flexible part / length 82.5 cm
280 174	ring-drainage-catheter Fr. 7.0 / 6 sideports / length 50 cm

disposable PTC-Set Fr. 8.4 / consisting of 7 parts

chiba needle



dilator




PTC needle




ring-drainage-catheter



guide wire



PTC exchange guide wire



drainage bag connector



Order No.	name
280 180	PTCD-Set Fr. 8.4 / consisting of 7 parts

single components:

280 181	chiba needle 22 G / length 20 cm
280 182	PTC needle 19 G incl. catheter Fr. 5.0 / length 25 cm
280 183	teflon coated guide wire .035 inch / 3 mm J-curved / length 80 cm
280 184	PTC exchanging guide wire .035 inch / 7.5 cm flexible part / length 82.5 cm
280 185	dilator Fr. 8.0 / length 20 cm
280 186	ring-drainage-catheter Fr. 8.4 / 32 sideports / length 50 cm
280 187	drainage bag connector / funnel-male / length 30 cm

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de