


STENT PLACEMENT SETS

re-usable



Stent Introduction Set **GPPS**

for placing of **Fr. 8.5 / Fr. 10.0 and Fr. 11.5** Stents in one step
(not to use for Pigtail and Double Pigtail Stents)



Advantage: Contrast injection while guide wire is inserted!



Order No.	diameter of external tube	fitting to stents	length	colour of external tube
200 080011	Fr. 8.5	Fr. 8.5	220 cm	orange
200 100011	Fr. 10.0	Fr. 10.0	220 cm	red
200 110011	Fr. 11.5	Fr. 11.5	220 cm	blue

Alternatively to GPPS: Pusher and Guiding Catheter separately

Pusher with distal radiopaque marker



Order No.	diameter	fitting to stents	length	colour of tube
200 05001	Fr. 5.0	Fr. 5.0	180 cm	dark grey
200 07001	Fr. 7.0	Fr. 7.0	180 cm	yellow
200 08001	Fr. 8.5	Fr. 8.5	180 cm	orange
200 10001	Fr. 10.0	Fr. 10.0	180 cm	red
200 11001	Fr. 11.5	Fr. 11.5	180 cm	blue

Guiding Catheter with distal radiopaque marker and detachable LL-connector



Order No.	diameter	fitting to stents	length	colour of tube	colour of LL-connector
200 00325	Fr. 5.0	Fr. 8.5	320 cm	white	blue
200 00326	Fr. 6.0	Fr. 10.0 / Fr. 11.5	320 cm	white	light grey

Legend:



This instrument is sterile packed separately.



This instrument is autoclavable by 134°C.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de