

EXTRACTION BALLOONS


disposable



 **ENDO-FLEX**[®]


disposable Extraction Balloons

- total length 200 cm
- 4 balloon sizes
- contrast injection port
- mandrin




double lumen

Order No.	Ø balloon deflated (mm)	Ø balloon inflated (mm)	air capacity (ml)	catheter	fitting to guide wire (inch)
134 000	1.9	9.0	1.5	Fr. 5.0	.025
134 005	1.9	12.0	2.5	Fr. 5.0	.025
134 010	2.6	14.0	2.5	Fr. 7.0	.035
134 011	1.9	14.0	2.5	Fr. 5.0 - Fr. 7.0	.035
134 020	2.6	18.0	4.0	Fr. 7.0	.035



triple lumen

Order No.	Ø balloon deflated (mm)	Ø balloon inflated (mm)	air capacity (ml)	catheter	fitting to guide wire (inch)
134 030	1.9	16.0	5.0	Fr. 7.0	.035
134 031	1.9	14.0	4.0	Fr. 7.0	.035
134 032	1.9	12.0	2.4	Fr. 7.0	.035
134 033	1.9	9.0	1.7	Fr. 7.0	.035



INSTRUMENT CONTAINS LATEX!

The latex and rubber devices should not be used in sensitive patients with allergic reaction to these substances. In these cases other alternative treatments are to be used.

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de