


kink resistant Memory Basket \varnothing 2.3 mm,
for working channels $\geq \varnothing$ 2.8 mm

basket made of gold plated Nitinol,
radiopaque tube,
incl. handle:




4 wires



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	146 224	146 225	146 226	141 224	141 225	141 226

6-fach



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	147 224	147 225	147 226	142 224	142 225	142 226

NOTE: Kink resistant Memory Baskets are not suitable for lithotripsy!

kink resistant Memory Basket \varnothing 1.8 mm,
for working channels $\geq \varnothing$ 2.0 mm

basket made of gold plated Nitinol,
radiopaque tube,
incl. handle:




4 wires



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
Nitinol wire	146 183	146 184	146 185	141 183	141 184	141 185

6 wires



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
Nitinol wire	147 183	147 184	147 185	142 183	142 184	142 185

NOTE: Kink resistant Memory Baskets are not suitable for lithotripsy!

Legend:

- This instrument is autoclavable by 134°C / 273°F.
- This instrument is rotatable.
- This instrument is detachable for cleaning.
- This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de

EXTRACTOR BASKETS
re-usable



© ENDO-FLEX - technical specification are subject to change

Extractor Baskets \varnothing 2.3 mm for working channels $\geq \varnothing$ 2.8 mm

basket made of stainless steel,
incl. handle:



4 wires




	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224	150 225	150 226	150 227	151 224	151 225	151 226	151 227
braided wire	160 224	160 225	160 226	160 227	161 224	161 225	161 226	161 227

6 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	152 224	152 225	152 226	152 227	153 224	153 225	153 226	153 227
braided wire	162 224	162 225	162 226	162 227	163 224	163 225	163 226	163 227

8 wires




	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	154 224	154 225	154 226	154 227	155 224	155 225	155 226	155 227
braided wire	164 224	164 225	164 227	164 227	165 224	165 225	165 226	165 227

Extractor Baskets \varnothing 2.3 mm, accepts guide wire .035 inch, for working channels $\geq \varnothing$ 2.8 mm

basket made of stainless steel,
incl. handle:




4 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224GW	150 225GW	150 226GW	150 227GW	151 224GW	151 225GW	151 226GW	151 227GW
braided wire	160 224GW	160 225GW	160 226GW	160 227GW				

6 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	152 224GW	152 225GW	152 226GW	152 227GW	153 224GW	153 225GW	153 226GW	153 227GW

recommenden equipment

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400




Extractor Baskets \varnothing 1.8 mm for working channels $\geq \varnothing$ 2.0 mm

basket made of stainless steel,
incl. handle:




4 wires



	dormia-type			basket-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
monofilament wire	150 183	150 184	150 185	151 183	151 184	151 185
braided wire	160 183	160 184	160 185	161 183	161 184	161 185

6 wires



	dormia-type			basket-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
monofilament wire	152 183	152 184	152 185	153 183	153 184	153 185