



Lithotripter Baskets \varnothing 2.6 mm for working channel $\geq \varnothing$ 2.8 mm

- Teflon[®] tube with lateral contrast injection port
- braided traction cable
- total length 4 m





4 wires



	basket type		
	basket size (in stretched condition)		
	50 mm	60 mm	70 mm
braided wire	122 265N	122 266N	122 267N

6 wires



	dormia type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	50 mm	60 mm	70 mm	50 mm	60 mm	70 mm
braided wire	124 265N	124 266N	124 267N	128 265N	128 266N	128 267N

recommended equipment

	handle	stainless steel spiral
Order No.	120 000N	for working channels $\geq \varnothing$ 2.8 mm: 120 261 (\varnothing 2.6 mm / length 195 cm)
Order No.	120 000N	for working channels $\geq \varnothing$ 3.7 mm: 120 331 (\varnothing 3.4 mm / length 195 cm)
Order No. (methodical)	120 500	for working channels $\geq \varnothing$ 2.8 mm: 120 527 (\varnothing 2.6 mm / length 210 cm)
Order No. (emergency)	120 500	without endoscope: 120 641 (\varnothing 4.0 mm / length 100 cm)

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de

MECHANICAL LITHOTRIPTORS

re-usable



Mechanicals Lithotripter Set Type I Ø 2.6 mm
handle with coaxial fixation, stainless steel spiral, 2 baskets


For methodical lithotripsy „Through-The-Scope“.



Teflon® tube with lateral contrast injection port



baskets made of braided wire, with 4 m traction cable



stainless steel spiral Ø 2.6 mm



handle with coaxial fixation

Set Type I	basket size (in stretched condition)		
	50 mm	60 mm	70 mm
Order No.	120 265N	120 266N	120 267N

Replacement

	handle	stainless steel spiral
Order No.	120 000N	for working channels ≥ Ø 2.8 mm: 120 261 (Ø 2.6 mm / length 195 cm)

Mechanicals Lithotripter Set Type II Ø 2.6 mm
handle and stainless steel spiral

Applicable with various baskets, suitable for methodical as well as emergency lithotripsies.



stainless steel spiral Ø 2.6 or 4.0 mm



handle

for methodical lithotripsy „Through-The-Scope“, stainless steel spiral 2.6 mm Ø, length 210 cm	
Order No.	120 526

for „Emergency“ lithotripsy without endoscope, stainless steel spiral 4.0 mm Ø, length 100 cm	
Order No.	120 640

Replacement

lithotripsy	handle	stainless steel spiral
Order No. (methodical)	120 500	for working channels ≥ Ø 2.8 mm: 120 527 (Ø 2.6 mm / length 210 cm)
Order No. (emergency)	120 500	without endoscope: 120 641 (Ø 4.0 mm / length 100 cm)