


## Achalasia Balloon Dilator

- total length 100 cm



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 3080	30.0 mm	80.0 mm		4.7 mm	.035 inch
341 3580	35.0 mm	80.0 mm		4.7 mm	.035 inch
341 4080	40.0 mm	80.0 mm		4.7 mm	.035 inch

### recommended equipment

Inflation / Deflation Device incl. PSI-Manometer 0-440 PSI



Order No. 340 800 INFL



kink resistant Nitinol Guide Wire „stripe guide“






length	.025 inch	.035 inch
400 cm	215 25400	215 35400



© ENDO-FLEX - Technical specification are subject to change

#### Legend:

-  This instrument is designed and warranted for single use only.
-  This instrument is sterile packed separately.
-  This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH  
 Alte Hünxer Str. 113 - 115  
 46562 Voerde / Germany  
 Phone: +49 281 - 9 44 00 - 0  
 Fax: +49 281 - 9 44 00 - 11  
 e-mail: export@endo-flex.de  
 www.endo-flex.de

# DILATION BALLOONS


disposable



 ENDO-FLEX®

## Biliary Balloon Dilator


- total length 200 cm
- radiopaque markers



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 0420	4.0 mm	20.0 mm	2.8 mm	2.3 mm	.035 inch
341 0620	6.0 mm	20.0 mm	2.8 mm	2.3 mm	.035 inch
341 0830	8.0 mm	30.0 mm	2.8 mm	2.3 mm	.035 inch

## Esophageal Balloon Dilator

- total length 200 cm
- radiopaque markers



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 0880	8.0 mm	80.0 mm	2.8 mm	1.9 mm	.035 inch
341 1280	12.0 mm	80.0 mm	2.8 mm	1.9 mm	.035 inch
341 1680	16.0 mm	80.0 mm	3.2 mm	1.9 mm	.035 inch
341 1880	18.0 mm	80.0 mm	3.2 mm	1.9 mm	.035 inch
340 0880	8.0 mm	80.0 mm	2.8 mm	1.9 mm	
340 1280	12.0 mm	80.0 mm	2.8 mm	1.9 mm	
340 1680	16.0 mm	80.0 mm	3.2 mm	1.9 mm	
340 1880	18.0 mm	80.0 mm	3.2 mm	1.9 mm	

### recommended equipment

Inflation / Deflation Device incl. PSI-Manometer 0-440 PSI



Order No. 340 800 INFL

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400

### recommended equipment

Inflation / Deflation Device incl. PSI-Manometer 0-440 PSI



Order No. 340 800 INFL

kink resistant Nitinol Guide Wire „stripe guide“



length	.025 inch	.035 inch
400 cm	215 25400	215 35400