

SPRAY CATHETERS

re-usable



radial Spray Catheter Ø 2.3 mm, for working channels ≥ Ø 2.8 mm Ø

- for dyeing suspicious mucosa by spraying colour
- 360° radial spraying
- LL connector with mandrin



length	Order No.
120 cm	47 226-A
160 cm	47 226-B
180 cm	47 226-C
230 cm	47 226-G

Spray Catheter for wide spray Ø 1.8 mm, for working channels ≥ Ø 2.0 mm

- for dyeing suspicious mucosa by spraying colour
- for straight, cone-shaped spraying
- LL connector with mandrin



length	Order No.
120 cm	47 185-A
160 cm	47 185-B
180 cm	47 185-C
230 cm	47 185-G

Legend:



This instrument is autoclavable by 134°C / 273°F.



This Instrument possesses a flushing port.

ENDO-FLEX GmbH
 Alte Hünxer Str. 113 - 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de