


# CLIP APPLICATION DEVICE

re-usable



# Clip Application Device $\varnothing$ 2.5 mm for working channels $\geq \varnothing$ 2.8 mm for placement of hemostasis and marking clips

- metal spiral in TEFLON<sup>®</sup> tube
- gold plated hook
- metal handle



length	Order No.
160 cm	D 2025-B
180 cm	D 2025-C
230 cm	D 2025-G

## Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is rotatable.



This instrument is detachable for cleaning.

ENDO-FLEX GmbH  
Alte Hünxer Str. 113 - 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)