

LOOP APPLICATION DEVICE

re-usable

AUTOCLEAVABLE
134°C
273°F



 **ENDO-FLEX**[®]

Loop Application Device \varnothing 2.5 mm for working channels $\geq \varnothing$ 2.8 mm for placement of Endo-Loops

- metal spiral in TEFLON[®] tube
- stationary mounted handle



length	Order No.
160 cm	2125-B
180 cm	2125-C
230 cm	2125-G

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de