

HF CURRENT DEVICES

re-usable

AUTOCLEARABLE
134°C
273°F



 **ENDO-FLEX[®]**

Handles with HF current connector re-usable

Universal Handle




Order-No.

62 0300



Handle for large Polypectomy Snare




Order No.

62 1300



Active Cords re-usable



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is designed to be used with HF current.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0

Fax: +49 281 - 9 44 00 - 11

e-mail: export@endo-flex.de

www.endo-flex.de