

COLON DECOMPRESSION PROBES

disposable



disposable Colon Decompression Probes

- accepts guide wire \leq .052 inch
- 6 distal sideports
- universal connector
- radiopaque tube



| tube | length | Order No. |
|----------|--------|-----------|
| Fr. 14.0 | 175 cm | 291 140 |
| Fr. 16.0 | 175 cm | 291 160 |
| Fr. 18.0 | 175 cm | 291 180 |

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de