


SUCTION POLYP TRAP™

disposable



 **ENDO-FLEX**®

Suction Polyp Trap™

for collecting polyps after excision

The Suction Polyp Trap™ processes:

- four filtered traps to collect and store polyps in different sections
- four open areas for straight unfiltered suction.

Two connection hoses at the lid:

- hose with connector to be attached to the suction device
- hose without connector to be attached to the endoscope.

The Polyp Trap™ can be sealed for transportation to a laboratory.



Order No.	package
900 334H	6 pieces

Legend:



This instrument is designed and warranted for single use only.

ENDO-FLEX GmbH
Alte Hünxer Str. 113 - 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de