




Useful Products

Order No.	Name
ZB 0030	E.R.C.P. cover cloth (1 package = 10 pieces)

- fixed to the trolley by adhesive stripes
- including 4 clear plastic folders for safe storage of devices and equipment



Legend:

-  This instrument is designed and warranted for single use only.
-  This instrument is sterile packed separately.
-  This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
 Alte Hünxer Str. 113 - 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de

© ENDO-FLEX - Technical specification are subject to change


USEFUL PRODUCTS



 **ENDO-FLEX**®

Useful Products

Order No.	Name
131 005	Biopsy Pincers




Order No.	Name
900 1514	Aspirator 25 ml (1 package = 10 pieces)



- incl. sealing lid and label

Order No.	Name
900 1515	Aspiration and Lavagesystem 200 ml (1 package = 20 pieces)



- tube lengths 150 cm
- Luer Lock connections for uncoupling and locking

Order No.	Name
199 551	Aqua-Gel



- non sterile lubricant on aqueous base
- anallergenic

Useful Products

Order No.	Name
199 520	silicone oil 50 ml
900 199 520	silicone oil 250 ml



- uncoloured silicone oil for valve care

Order No.	Name
ZB 0010	Colonoscopy Shorts (1 package = 10 pieces)



- three-ply material, water-repellent, bacteria-restraining, breathe-active
- Universal size
- slit on back side

Order No.	Name
ZB 0018	Gastroskopy Napkin (1 package = 10 pieces)



- safe absorption of saliva and secretion