

ENDOCART STANDARD ES-100

- Durable plastic work surface
- Two adjustable shelves
- 6" medical casters (two locking)
- Fold-down keyboard tray with mouse pad
- Hospital grade power strip
- Two stainless steel handles
- Easy access rear cover
- Side panel colour options (red, green, blue)
- Ready to receive optional monitor arm mounts, locking drawer, dual scope hanger and more.



ENDOCART BASIC EB-100

- Durable plastic work surface
- Adjustable shelf
- 5" medical casters (two locking)
- Keyboard tray
- Basic power strip
- Ready to receive flat panel arm mounts
- Ready to receive optional locking drawer



ENDOCART VARIABLE HEIGHT ROLL STAND ERS-100

- Adjustable height (47-62.5 inches)
- Standard flat panel arm
- 4" casters (two locking)
- 21 LB weighted base
- Optional wire basket



ENDOCART MONITOR STAND EMS-104

- 4" medical casters (two locking)
- Stainless steel handle
- Base shelf platform
- Holds flat panel and CRT monitors
- Ready to dock into EndoCart Deluxe
- Ready to receive optional shelves (EMS-201)
- Ready to receive optional drawers

