


ENDOCART® DELUXE ED-100

- Large durable plastic work surface
- Three drawers (1 locking) or Two drawers (1 large)
- Two adjustable shelves
- 6" medical casters (2 locking, 2 directional)
- Fold-down keyboard tray with mouse pad
- Hospital grade power strip
- 2 stainless steel handles
- Easy access rear covers
- Side panel color options (red, green, blue)
- Ready to receive docking Monitor Stand
- Ready to receive optional monitor arm mounts that can be positioned anywhere on the back of the cart

Color Panel Choices:



The Ideal Mobile System



The EndoCart Deluxe and the Monitor Stand easily lock together to move as one system!